



# PESI EREE2

## WEIGHT SET WITH WOODEN CASE AND STAINLESS STEEL WEIGHTS IN E2 CLASS



Weight set with wooden case, fitted with stainless steel weights in class E2 according to OIML recommendation R111 and suitable for legal metrology and industrial applications. The precision class E2 allows the use of these weights for checking the calibration of scales and instruments in class I up to 500.000e.

Ideal for the verification of weighing instruments in ISO quality regime.


Our solutions for quality assurance and certainty of your measurements.

## MAIN FEATURES WEIGHTS

- Polished and austenitic stainless steel weights
- E2 precision class
- Compliant with OIML recommendation - R111
- Density 7950 kg/m<sup>3</sup>
- Accredia (EA) calibration certificate upon request
- Initial verification stamp available upon request

## MAIN FEATURES CASE

- Wooden
- Fitted with glove and/or tweezer (tweezer and dust brush with the WSE2M05)



## LASER MARKING SERVICE UPON REQUEST:

- The laser marking service is useful to identify each weight or mass individually. More specifically:
  - The marking allows easier identification and traceability of each weight or mass;
  - It is perfectly smooth, does not remove material from the weight and does not affect the surface, thus avoiding dust deposit;
  - It is customisable, according to the requirements of the companies.

In addition:

- It remains unchanged over time;
- It can be performed on previously purchased masses and weights, of all accuracy classes, from 1 g to 20 kg;
- It is perfect for pharmaceutical industries and laboratories.

### CALIBRATION SERVICE

Below you can see a preview of the calibration certificate. Please pay special attention to point 2. "Required information", data to be communicated to the laboratory when requesting calibration.

Calibration Certificate

## DETAIL 1



Example of optional laser marking

## WEIGHT SETS E2 CLASS

Codice	Numero pezzi	Pesi contenuti	Note
WSE2M05	12	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	Total: 1,11g
WSE2M50	20	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50	Total: 111,11g
WSE2M100	21	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100	Total: 211,11g
WSE2M200	23	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200	Total: 611,11g
WSE2M500	24	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	Total: 1111,11g
WSE2MK1	25	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000	Total: 2111,11g
WSE2MK2	27	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000	Total: 6111,11g
WSE2MK5	28	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000	Total: 11111,11g
WSE2MK10	29	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000, 10000	Total: 21111,11g
WSE250	8	g: 1, 2, 2, 5, 10, 20, 20, 50	Total: 110g
WSE2100	9	g: 1, 2, 2, 5, 10, 20, 20, 50, 100	Total: 210g
WSE2200	11	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200	Total: 610g
WSE2500	12	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	Total: 1110g
WSE2K1	13	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000	Total: 2110g
WSE2K2	15	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000	Total: 6110g
WSE2K5	16	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000	Total: 11110g
WSE2K10	17	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000, 10000	Total: 21110g



QUALITY ASSURANCE AND CERTAINTY OF YOUR MEASUREMENTS



**OFFICE**

Via Picasso, 18/20  
20025 Legnano (MI) - Italy

SALES AND TECHNICAL ASSISTANCE SERVICE

*All the specifications and images can be subject to variations and upgrades without prior notice.*