



## MASSE A DISCO

### M1 PRECISION CLASS DISK WEIGHTS AND HOLDING BARS



Disc weights suitable for legal metrology applications and the industrial sector and fully customisable. They have a groove and relief to allow them to interlock and keep them concentric; the slot allows them to be positioned on the bar. Their tolerances are those provided for OIML M1 class masses and therefore allow these weights to be used for the control and calibration of scales and class III instruments up to 10.000e. Ideal for ISO quality weighing instrument verification.

Our solutions for quality assurance and certainty of your measurements.

## MAIN FEATURES

- M1 PRECISION CLASS
- MADE IN PHOSPHATE IRON
- ACCREDIA (EA) CALIBRATION CERTIFICATE, UPON REQUEST
- ALL MASSES ARE FITTED WITH A CAVITY AND RELIEF TO ALLOW THEIR INTERLOCKING AND KEEP THEM CONCENTRIC
- THE DISC MASSES ARE FITTED WITH A SLOT FOR THE POSITIONING ON THE BARS
- THE DIMENSIONS INDICATES ARE CUSTOMISABLE.

### CALIBRATION SERVICE

Below you can see a preview of the calibration certificate. Please pay special attention to point 2. "Required information", data to be communicated to the laboratory when requesting calibration.

Calibration Certificate



QUALITY ASSURANCE AND CERTAINTY OF YOUR MEASUREMENTS



**OFFICE**

Via Picasso, 18/20  
20025 Legnano (MI) - Italy

SALES AND TECHNICAL ASSISTANCE SERVICE

*All the specifications and images can be subject to variations and upgrades without prior notice.*