



PIPETTE

CALIBRATION OF OTHER INSTRUMENTS: PIPETTES



ACCREDIA calibration service of pipettes (microdosers).

Our solutions for quality assurance and certainty of your measurements.

PIPETTES

- ACCREDIA calibration of fixed or variable volume, single- or multi-channel pipettes in accordance with ISO 8655:2022.
- Measurements are taken on 10 dispensing (fixed volume) or 10 dispensing on 3 points (variable volume).
- The LAT calibration certificate is issued.
- The service includes pipette maintenance carried out before the calibration.
- Preventive maintenance consists of greasing and replacing wear parts with original spare parts.
- Calibration service time: 8 days from receipt of order.
- Service can also be provided at the customer's site (non-accredited activity).

- CALIBRATION SERVICE

Below you can see a preview of the calibration certificate. Please pay special attention to point 2. "Required information", data to be communicated to the laboratory when requesting calibration.

Calibration Certificate



QUALITY ASSURANCE AND CERTAINTY OF YOUR MEASUREMENTS



OFFICE

Via Picasso, 18/20
20025 Legnano (MI) - Italy

SALES AND TECHNICAL ASSISTANCE SERVICE

All the specifications and images can be subject to variations and upgrades without prior notice.