





PESO_MINIMO



MINIMUM WEIGHT CERTIFICATE



Applicable for all types of scales.

Our solutions for quality
assurance and certainty
of your measurements.

- 
- 
- When choosing the right balance for weighing a specific sample, resolution and linearity data alone can be limiting. They, in fact, do not consider the uncertainty of use of the balance.

- 
- 
- So, it is necessary to include another aspect, which is the minimum weight, determined according to considerations inside the Euromet cg-18 guide and that represents the smallest quantity of a sample that is possible to weigh on a specific scale obtaining a reliable result.

- 
- 
- Download the brochure for more information:

[Minimum weight certificate](#)

Request now the minimum weight certificate for your scale!



QUALITY ASSURANCE AND CERTAINTY OF YOUR MEASUREMENTS



OFFICE

Via Picasso, 18/20
20025 Legnano (MI) - Italy

SALES AND TECHNICAL ASSISTANCE SERVICE

All the specifications and images can be subject to variations and upgrades without prior notice.