



PESI EREM1\_I NOX

WEIGHT SET WITH WOODEN CASE AND STAINLESS STEEL WEIGHTS IN M1 CLASS



Weight set with wooden case, fitted with stainless steel weights in M1 class according to OIML recommendation R111 and suitable for legal metrology and industrial applications. The precision class M1 allows the use of these weights for checking and calibrating scales and instruments in class III up to 10.000e.

Ideal for the verification of weighing instruments in ISO quality regime.


Our solutions for quality assurance and certainty of your measurements.

## MAIN FEATURES WEIGHTS

- Polished and austenitic stainless steel weights
- M1 precision class
- Compliant with OIML recommendation - R111
- Density 7950 kg/m<sup>3</sup>
- Accredia (EA) calibration certificate upon request
- Initial verification stamp available upon request

## MAIN FEATURES CASE

- Wooden
- Complete with glove and/or tweezers and/or brush



## LASER MARKING SERVICE UPON REQUEST:

- The laser marking service is useful to identify each weight or mass individually. More specifically:
    - The marking allows easier identification and traceability of each weight or mass;
    - It is perfectly smooth, does not remove material from the weight and does not affect the surface, thus avoiding dust deposit;
    - It is customisable, according to the requirements of the companies.
- In addition:
- It remains unchanged over time;
  - It can be performed on previously purchased masses and weights, of all accuracy classes, from 10 mg to 20 kg;
  - It is perfect for pharmaceutical industries and laboratories.

### CALIBRATION SERVICE

Below you can see a preview of the calibration certificate. Please pay special attention to point 2. "Required information", data to be communicated to the laboratory when requesting calibration.

Calibration Certificate


## DETAIL 1



Example of optional laser marking

## VERSIONS

### WEIGHT SETS M1 CLASS

Code	Number of pieces	Weights inside	Notes
 WSAM1M05	12	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	Total: 1,11g Complete with tweezers and brush
WSAM1M50	20	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50	Total: 111,11g
 WSAM1M100	21	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100	Total: 211,11g
WSAM1M200	23	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200	Total: 611,11g
WSAM1M500	24	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	Total: 1111,11g
WSAM1MK1	25	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000	Total: 2111,11g
 WSAM1MK2	27	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000	Total: 6111,11g
 WSAM1MK5	28	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000	Total: 11111,11g
 WSAM1MK10	29	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000, 10000	Total: 21111,11g
 WSAM150	8	g: 1, 2, 2, 5, 10, 20, 20, 50	Total: 110g
 WSAM1100	9	g: 1, 2, 2, 5, 10, 20, 20, 50, 100	Total: 210g
WSAM1200	11	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200	Total: 610g
WSAM1500	12	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	Total: 1110g
WSAM1K1	13	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000	Total: 2110g



WSAM1K2

15

g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000

Total: 6110g



WSAM1K5

16

g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000

Total: 11110g



WSAM1K10

17

g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000,  
5000, 10000

Total: 21110g



QUALITY ASSURANCE AND CERTAINTY OF YOUR MEASUREMENTS



**OFFICE**

Via Picasso, 18/20  
20025 Legnano (MI) - Italy



SALES AND TECHNICAL ASSISTANCE SERVICE

*All the specifications and images can be subject to variations and upgrades without prior notice.*