



PESI EREE2\_PLASTICA

WEIGHT SET WITH PLASTIC CASE AND STAINLESS STEEL WEIGHTS IN E2 CLASS



Weight set with high-quality plastic case, complete with stainless steel weights in class E2 according to OIML R111 recommendation and suitable for legal metrology and industrial applications. The E2 accuracy class allows these weights to be used for checking and calibrating class I scales and measuring instruments up to 500,000e.

Ideal for checking weighing instruments in accordance with ISO quality standards.


Our solutions for quality assurance and certainty of your measurements.

## MAIN FEATURES WEIGHTS

- Polished and austenitic stainless steel weights
- E2 precision class
- Compliant with OIML recommendation R111
- Density 7950 kg/m<sup>3</sup>
- Accredia (EA) calibration certificate available on request
- Verification stamp available on request

## MAIN FEATURES CASE

- Made of high-quality plastic
- Complete with glove and/or tweezers and/or brush



## LASER MARKING SERVICE UPON REQUEST:

- The laser marking service is useful to identify each weight or mass individually. More specifically:
    - The marking allows easier identification and traceability of each weight or mass;
    - It is perfectly smooth, does not remove material from the weight and does not affect the surface, thus avoiding dust deposit;
    - It is customisable, according to the requirements of the companies.
- In addition:
- It remains unchanged over time;
  - It can be performed on previously purchased masses and weights, of all accuracy classes, from 10 mg to 20 kg;
  - It is perfect for pharmaceutical industries and laboratories.

### CALIBRATION SERVICE

Below you can see a preview of the calibration certificate. Please pay special attention to point 2. "Required information", data to be communicated to the laboratory when requesting calibration.

Calibration Certificate

DETAIL 1



## VERSIONS

### E2 CLASS PLASTIC WEIGHT SET

Code	Number of pieces	Weights inside	Notes
 WSE2M50PL	20	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50	Total: 111,11g
 WSE2M100PL	21	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100	Total: 211,11g
 WSE2M200PL	23	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200	Total: 611,11g
 WSE2M500PL	24	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	Total: 1111,11g
 WSE2MK1PL	25	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000	Total: 2111,11g
 WSE2MK2PL	27	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000	Total: 6111,11g
 WSE2MK5PL	28	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000	Total: 11111,11g
 WSE250PL	8	g: 1, 2, 2, 5, 10, 20, 20, 50	Total: 110g
 WSE2100PL	9	g: 1, 2, 2, 5, 10, 20, 20, 50, 100	Total: 210g
 WSE2200PL	11	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200	Total: 610g



WSE2500PL

12

g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500

Total: 1110g



WSE2K1PL

13

g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000

Total: 2110g



WSE2K2PL

15

g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000

Total: 6110g



WSE2K5PL

16

g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000

Total: 11110g



QUALITY ASSURANCE AND CERTAINTY OF YOUR MEASUREMENTS



**OFFICE**

Via Picasso, 18/20  
20025 Legnano (MI) - Italy



SALES AND TECHNICAL ASSISTANCE SERVICE

*All the specifications and images can be subject to variations and upgrades without prior notice.*