



CAMPIONATORI ARIA

OTHER INSTRUMENTS CALIBRATION: AIR



Calibration service of microbiological air samplers, performed with LAT certified instrumentation (or equivalent in EA). Ideal for instrumentation used in the pharmaceutical, biotechnology and food industries.

Our solutions for quality assurance and certainty of your measurements.

AIR SAMPLERS

- Calibration of microbiological air samplers.
- Measurement of air velocity (m/s) during a 500 l cycle ("As found" calibration). Further 500 l cycle and recording of the time taken to complete it ("As left" calibration).

Calibration service time: 10 days from receipt of order.



QUALITY ASSURANCE AND CERTAINTY OF YOUR MEASUREMENTS



OFFICE

Via Picasso, 18/20
20025 Legnano (MI) - Italy



SALES AND TECHNICAL ASSISTANCE SERVICE

All the specifications and images can be subject to variations and upgrades without prior notice.